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Operating Manual Conductivity Meter LWM S3

Werner Order No.: 1712-046

Description:

- 1 Display (digital)
- 2 Rotary switch for switch point
- 3 LED (red / green)
- 4 Connection electrode
- 5 Connection power supply

The conductivity meter LWM S3 provides a permanent observation of the pure water quality of ion exchanger columns. The measuring value is digitally displayed in $\mu\text{S}/\text{cm}$, the unit of electrical conductivity. A LED, placed on the front plate, shows green light, as long as the current limit value complies. If the current limit value is exceeded, the LED turns to red light. If the measuring range is exceeded, "OFL" will be displayed.

Initiation:

4. If necessary, place electrode on ion exchanger column before. See user instruction / set into operation for AQUADEM® DL, -DF respectively -SDF.
5. Plug conductivity meter on the electrode connection on the ion exchanger column.
6. Connect power supply to conductivity meter and plug into a socket.

Set switch point:

For switch point setting, a rotary switch is placed at the left side of the box. By using a screwdriver, the required switch point can be set. During the setting procedure the selected switch point value is displayed for 3 s.

Technical Data:

Measuring range	Voltage power supply	Voltage device	Conductivity cell factor	Switch points
0-199,9 $\mu\text{S}/\text{cm}$	100-240V AC	9V DC	$C = 0,1 \text{ cm}^{-1}$	0,5 / 1,3 / 2 / 4,3 / 7,5 / 10 / 15 / 20 / 30 / 50 $\mu\text{S}/\text{cm}$